

## **DISINFECTING CLEANER**

Cleaner & Disinfectant for Spray and Wipe Cleaning in Schools, Health Care and Other Facilities

Generate for 15 minutes; Disinfect for 14 days 2.4oz of Saline + 0.5 Gallon H<sub>2</sub>0 disinfects approx. 4000<sup>2</sup>ft

The active ingredients or mixed oxidants generated in GenEon Solution have been shown to be effective against a broad range of blood borne pathogens including:

- **Hepatitis A**
- **Hepatitis C**
- **Hepatitis B**
- HIV
- **Herpes Simplex Virus Type 1 & 2**

Effective when used as directed on hard non-porous inanimate surfaces such as: Toilet Bowl, Restrooms, Chairs, Kitchen Areas, Chopping Boards, Tables, Blood & Organic Matter, Door Knobs, Counter Tops, Stainless Steel, Sinks, Natural Stone, and other similar surfaces. \*Test in an inconspicuous area before using.

The same			에 트로피스 파스 등 글로 사내가 되었다.
-	pH 5.5-7.5 FAC ppm	30 Seconds	
	<199	Listeria Monocytogenes E.coli O157.H7 Salmonella enterica	
	425	Bacteriophage MS-2 Pseudomonas aeruginosa Mycobacterium terrea (TB)	EFFICACY 1
96	600	Campylobacter Enterovirus+ Legionella Hepatitis B	Kills <u>Bacteria</u> , <u>Vir</u> In as Little as 30  FAC= Free Available Chlo MRSA = 2.5 Minutes Conf
	950-1100	HIV	
	1100-1200	Murine Noro Virus C. diffcile++ Polio VirusCHAT Lsc1 Hepatitis C Herpes	



GenEón

